

Genotyping of Mouse Tail DNA via PCR

I. *Mouse tailing*

[Pups are tailed (for DNA) and toed (for identification) between 8-14 days of age.]

- A. Remove tail sample of approximately 0.25 inches by pinching the tip of the tail to expel blood and cutting with scissors.
- B. Place tail sample in 1.5 mcf tube for digestion.
- C. Scissors is washed liberally with 75% EtOH between cuts to prevent cross contamination.

II. *Tail lysis and sample preparation*

Preparation of tail lysis buffer

Final concentration [100mM Tris-HCl, pH8.8; 5mM EDTA, pH8.0; 0.2%SDS; 200mMNaCl]

- **Just prior to use, add proteinase K (Stock conc= 20mg/mL) to the buffer to achieve a final concentration of 100 ug/mL.**
- A. Add 500 uL of tail lysis buffer to each tube and vortex briefly to mix.
 - B. Tails are incubated o/n at 55C with agitation (250rpm on a heated shaker).
 - C. The following morning, tubes are vortexed and centrifuged [5 min x 14000 rpm@RT] to pellet debris.
 - D. Dilute lysis 1:11 in sterile H₂O for use in PCR analysis.

III. *PCR*

A. Primer sequences and product size

| GENE | Primer name | Primer sequence | Amplicon size |
|------------------|-------------|-------------------------|---------------|
| | | | |
| <i>nfi-a</i> | I2B | TGCTGTGTTCTGGTCAGTCAAG | 405 bp |
| | I2CC | CAAAGCAAATCTCCATGCTCGG | |
| | | | |
| <i>neo</i> | neo57 | GGAGAGGCTATTCGGCTATGAC | 315 bp |
| | neo371R | CGCATTGCATCAGCCATGATGG | |
| | | | |
| <i>sry</i> (sex) | sry1 | AACAAGTGGGCTTTGCACATTG | 166, 146 bp |
| | sry2 | GTTTATCAGGGTTTCTCTCTAGC | (doublet) |

B. PCR Reaction mixture

| | COMPONENT | | [STOCK] | [FINAL RXN] |
|-----------|--------------------------------------|-------------------|------------------|--------------------|
| 1. | 10X <i>Taq</i> Buffer | | 10X | 1X |
| | | Tris-HCl | 100mM | 10mM |
| | | MgCl ₂ | 15mM | 1.5mM |
| | | KCl | 500mM | 50mM |
| | | BSA | 0.2mg/mL | 0.02mg/mL |
| 2. | Multiplex primer mix | | 10X | 1X |
| | | I2B | 2 uM | 200 nM |
| | | I2CC | 2 uM | 200 nM |
| | | neo57 | 2 uM | 200 nM |
| | | neo371R | 2 uM | 200 nM |
| | | sry1 | 2 uM | 200 nM |
| | | sry2 | 2 uM | 200 nM |
| 3. | dNTP mix | | 1 mM each | 200 uM each |
| | (dA,dT,dG,dC) | | | |
| 4. | <i>taq</i> polymerase | | 5 U/uL | 0.5 Units |
| | (Life Technologies Cat# 10966026) | | | |
| | | | | |
| | | | | |
| | q.s. to final volume | | | 19.0 uL |
| | With H ₂ O | | | |
| | | | | |
| | DNA (diluted lysate) | | | 1.0 uL |
| | | | | |
| | | | | |
| | | | | 20.0 uL rxn |
| | | | | |

C. PCR Parameters

1. 94C for 4 min (initial melt)
2. 94C, 60C, 72C x 33 cycles (amplification)
3. 72C for 9 min (final extension)

D. PCR products are separated on a 2% agarose gel, stained with EtBR, and photographed for documentation.