

Induction of macrophage migration into peritoneal cavity

Prepare Thioglycollate Medium Brewer Modified (BBL, Becton Dickinson Microbiology Systems, MD, cat nr 211716): 38.5 g in 1L. Boil for 1 min to dissolve, and autoclave.

Inject the mice with 1 ml IP, and leave them for 3-5 days.

Collection of macrophages

Reagents:

1. Syringes with 5 ml sterile PBS and 18G needles.
2. Tubes with 2ml HBSS w/o Ca^{2+} or Mg^{2+} +10% FBS+2% Penicillin / Streptomycin, on ice.
3. Tube with Ketamine/Xylazine (K/X) at 10 mg/ml K, 2mg/ml X, in PBS.

Inject the mice with 100 μl K/X intramuscular (1 mg K, 0.2 mg X).

When KO, hold them with rubber bands and inject 5ml PBS in the belly, about 0.5 cm from sternon. Collect 3 ml from there (or 4 if the liquid is too clear). Place the liquid in the tubes without making bubbles.

Centrifuge tubes 5 min 1200 rpm.

Remove SN and add 1 ml DMEM+5%FBS+2% Pen/Strep (P 200U/ml, S 200 $\mu\text{g}/\text{ml}$).

Resuspend thoroughly and transfer 20 μl to 10ml counting medium.

Count and dilute to $1 \cdot 10^6$ cells/ml (If $N=6.2 \cdot 10^6$ then just need to add 5.2 ml DMEM; If $N=12.3$ ml, just add 11.3 ml DMEM, etc).

Place 50 μl of this dilution on columns 1 and 2 of a flat bottom (Costar) 96 well plate. The next mouse will use cols 3 and 4, next 5 and 6... ATTENTION: Cols 11 and 12 must go empty, since they are controls.

Incubate for at least 20 min at 37 °C.

Meanwhile, prepare the TLR3 agonists: Poly IC (5 $\mu\text{g}/\text{ml}$, Amersham)

After the incubation period for the M Φ is over, shake the plates to lift lymphocytes and remove the supernatant. Add 100 μl of the agonist solutions. Incubate for 4 h at 37°C.

Transfer the SN to a clean plate and freeze that at -80°C.

MTT marking of Macrophages

Immediately, add 50 μl of MTT (stock 2 mg/ml in PBS) diluted 1:4 in DMEM+10%FBS+2% Pen/Strep to the macrophages. Incubate from 1 h to O/N.

L929 cells

Culture medium: DMEM+10%FBS+2%Pen/Strep

To prepare bioassay plates, rinse an 80% confluent flask with PBS. Add 3 ml Trypsine EDTA and incubate until detachment. Add 8 ml medium, pool all flasks and centrifuge 5 min 1200 rpm.

Discard supernatant and add 10-20 ml medium.

Transfer 20 μ l of cells to 10 ml counting solution, and count.

TO PASS CELLS: transfer 200 μ l of the cells to a large flask.

TO PREPARE BIOASSAY: Dilute to $5 \cdot 10^5$ cells/ml, and pipette 100 μ l per well of a 96 well plate. Incubate 24 h.

Bioassay

Add 50 μ l Cyclohexemide (CHX) (From stock at 10 mg/ml or 10,000 U/ml, dilute 1:25 in DMEM+ 10%FBS+ 2%P+S) and Incubate 20 min 37 °C.

In the meantime, prepare TNF- α standard curve: Stock aliquoted at 10 μ l, corresponding to 100 ng. So, add 1 ml of medium to dilute it to 100ng/ml, from there, transfer 5 μ l to 1 ml of medium (0.5 ng/ml). Prepare standard on a 96 well plate: From dilution at 0.5 ng/ml, add 30 μ l to 270 μ l of medium (line A). From there down, make a serial dilution 150 μ l+150 μ l medium until line G. REMEMBER, line H is only medium!.

From the standard prepared, add 50 μ l on cols 11 and 12 (from line A to A 11 and A 12; from line B to B 11 and B 12, etc).

To the rest of the plate (cols 1-10), add 46 μ l of DMEM+10%FBS+2%P+S and 4 μ l of supernatant from M Φ .

Incubate the plates O/N at 37 °C.

Following day, add 70 μ l of pure MTT to the wells and incubate for at least 3h.

Dry and leave at RT.

To measure, lyse cells with 100 μ l lysis solution (2-propanol+10% H₂O+0.5% SDS+0.04N HCl: 3,6 L isopropanol+400 ml H₂O+20g SDS+4 ml HCl 4N) and incubate at RT for at least 30 min. Read at 590 nm using a reference filter of 690 nm.