




Welcome to The Genotyping Protocol System

Master Protocol


Strain Name: C57BL/6NJ-Prss36^{em1J}J
Stock Number: 029276
Allele: Prss36^{em1J}
Protocol Name: Prss36^{em1J}
Method: SEPARATED MELT
Created: 03-June -2016 (JKELMEN) **Updated:** 10-August -2016 (ESCHAAB)

Notes


Expected Results: Mutant = 232 bp
 Heterozygote = 140 bp and 232 bp
 Wild type = 140 bp

Attachments

View 29276 Prss36 gel.jpg



View 29276 Prss36 mut.jpg



View 29276 Prss36 wt.jpg

- View Prss36em1jMay2016.docx
- View Prss36 genomic.gcs

Protocol Primers

Primer	5' Label	Sequence 5' --> 3'	3' Label	Description	Reaction
28306	-	TCA TGA GTG AAT TGG AAG TGC	-	Common	A + B
28307	-	CTC CTC TCC TAA TCC CTT GGA	-	Wild type Reverse	A
28308	-	CGG GGA GCT AGG GTC AGT AG	-	Mutant Reverse	B

Number Of Reactions 1

Reaction A (2 Primer)		Cycling (TouchDown 65-60_12-12-11)			
Component	Final Concentration	Step #	Temp°C	Time	Note
ddH2O		1	94	2 min	-
Kapa 2G HS buffer	1.3 X	2	94	20sec	-
MgCl2	2.6 mM	3	65	15sec	-0.5 C per cycle decrease
dNTPS-kapa	.26 mM	4	68	10sec	-
28307	.5 uM	5	-	-	repeat steps 2-4 for 10 cycles (Touchdown)
		6	94	15sec	-

Reaction A (2 Primer)		Cycling (TouchDown 65-60_12-12-11)			
28306	.5 uM	7	60	15sec	-
Glycerol	6.5 %	8	72	10sec	-
Dye	1 X	9	-	-	repeat steps 6-8 for 28 cycles
Kapa 2G HS taq polymerase	.03 U/ul	10	72	2 min	-
DNA		11	10	-	hold

Reaction B (2 Primer)		Cycling (TouchDown 65-60_12-12-11)			
Component	Final Concentration	Step #	Temp°C	Time	Note
ddH2O		1	94	2 min	-
Kapa 2G HS buffer	1.3 X	2	94	20sec	-
MgCl2	2.6 mM	3	65	15sec	-0.5 C per cycle decrease
dNTPS-kapa	.26 mM	4	68	10sec	-
28306	.5 uM	5	-	-	repeat steps 2-4 for 10 cycles (Touchdown)
28308	.5 uM	6	94	15sec	-
Glycerol	6.5 %	7	60	15sec	-
Dye	1 X	8	72	10sec	-
Kapa 2G HS taq polymerase	.03 U/ul	9	-	-	repeat steps 6-8 for 28 cycles
DNA		10	72	2 min	-
		11	10	-	hold

Version 3.2