

**GENOTYPING BY PCR PROTOCOL
MUTANT MOUSE REGIONAL RESOURCE CENTER**

sacoord@mmrrc.org

800-910-2291 North America, +1-530-757-5710 International

Please provide the following information required for genetic analysis of your mutant mice.

Note to MAC users: to ensure your graphic can be viewed on a PC please follow the steps below when inserting the graphic into this document. DO NOT drag and drop or copy/paste the graphic into this document.

- Open the original graphic in the program that created it
- Choose File, Save As
- Select No Compression in the save options.
- Save as JPG or PNG or similar format that's compatible with both PC and Mac Word versions.
- Switch to Word, choose Insert, Picture, From File and choose the newly saved picture.

These instructions are very generic. The menu options for your graphics program may be different.

Donating Investigator/PI		
Dr. Mansoor Husain		
Email		
Mansoor.husain@uhn.ca		
Institution		
Toronto General Hospital Research Institute, University Health Network, Canada		
Address		
Room # PMCRT-3-910F, 101 College Street, MaRS East Tower,		
City	State	Zip
Toronto	ON	M5G 1L7
Lab Contact		
Dr. Talat Afroze, PhD		
Email		
Talat.afroze@utoronto.ca		
Telephone	FAX	
416-581-7491	NA	
Strain Name		MMRRC Stock Number
B6.Cg-Tg(tetO/CMV-ATP2B4,-lacZ)1Husa Tg(Myh6-tTA)6Smbf/Mmjax		MMRRC : 43951

GENOTYPING BY PCR PROTOCOL
MUTANT MOUSE REGIONAL RESOURCE CENTER

sacoord@mmrrc.org

800-910-2291 North America, +1-530-757-5710 International

NAME OF PCR: human PMCA4b Genomic PCR **MMRRC:** 043951

Protocol: *(PCR protocol provided by Donating Investigator)*

Reagent/Constituent	Volume (µL)
Water autoclaved; MilliQ brand reverse osmosis water	12.5
10x Buffer NH ₄ SO ₄ based Buffer; ThermoFisher	2.5
MgCl ₂ (stock concentration is 25mM) 2.5mM working concentration	2.0
Betaine (stock concentration is 5M) <i>Optional</i> N.A.	0
dNTPs (stock concentration is 10mM) 0.2mM working concentration	0.5
DMSO <i>Optional</i> N.A.	0
Human PMCA4b Forward Primer (stock concentration= 20µM)	0.5
Human PMCA4b Reverse Primer (stock concentration is 20µM)	0.5
	0
	0
Taq Polymerase 1Units/µL use 1.5Units per 25uL PCR reaction	1.5
DNA (5uL of crude tail digest using Proteinase K overnight; volume 400uL per tail)	5
<i>The total volume is auto-calculated based on volumes entered, right click the total and update field to show/recalculate the total volume.</i>	TOTAL VOLUME OF REACTION: 25.000 µL

Comments on protocol:

- We use ThermoFisher recombinant Taq DNA Polymerase 1U/uL Cat. # EP0404

Strategy: human PMCA4b Genomic PCR

Steps	Temp (°C)	Time (m:ss)	# of Cycles
1. Initiation/Melting HOT START? Yes X	94	3 min	1
2. Denaturation	94	30 sec	
3. Annealing steps 2-3-4 cycle in sequence	60	45 sec	45x
4. Elongation	72	60 sec	
5. Amplification	72	7 min	1
6. Finish	4	∞	n/a

Primers: huPMCA4b

Electrophoresis Protocol:

Name	Nucleotide Sequence (5' - 3')	Argarose: 2% in v:110		
1. huP4 For	GGCTCCCTGAGTGTACTCCC	Estimated Running:Time: 60 min.		
2. huP4 Rev	CCTGATGACGGTGCTCATTG	Primer Combination	Band	Genotype
3.		huP4 For + Rev	554 bp	huPMCA4b
4.			bp	
5.			bp	

NAME OF PCR: Tet O (Tet-sensitive TransActivator; tTA) Genomic PCR **MMRRC:** 043951

PCR protocol provided by Donating Investigator

**GENOTYPING BY PCR PROTOCOL
MUTANT MOUSE REGIONAL RESOURCE CENTER**

sacoord@mmrrc.org

800-910-2291 North America, +1-530-757-5710 International

Protocol: *(PCR protocol provided by Donating Investigator)*

Reagent/Constituent	Volume (μL)
Water autoclaved; MilliQ brand reverse osmosis water	12.5
10x Buffer NH ₄ SO ₄ based Buffer; ThermoFisher	2.5
MgCl ₂ (stock concentration is 25mM) 2.5mM working concentration	2.0
Betaine (stock concentration is 5M) <i>Optional</i> N.A.	0
dNTPs (stock concentration is 10mM) 0.2mM working concentration	0.5
DMSO <i>Optional</i> N.A.	0
Tet O (tTA) Forward Primer (stock concentration= 20μM)	0.5
Tet O (tTA) Reverse Primer (stock concentration is 20μM)	0.5
	0
	0
Taq Polymerase 1Units/μL use 1.5Units per 25uL PCR reaction	1.5
DNA (5uL of crude tail digest using Proteinase K overnight; volume 400uL per tail)	5
<i>The total volume is auto-calculated based on volumes entered, right click the total and update field to show/recalculate the total volume.</i>	TOTAL VOLUME OF REACTION: 25.000 μL

Comments on protocol:

- We use ThermoFisher recombinant Taq DNA Polymerase 1U/uL Cat. # EP0404

Strategy: Tet O (Tet-sensitive TransActivator; tTA) Genomic PCR

Steps	Temp (°C)	Time (m:ss)	# of Cycles
1. Initiation/Melting HOT START? Yes X	94	3 min	1
2. Denaturation	94	30 sec	
3. Annealing steps 2-3-4 cycle in sequence	50	45 sec	45x
4. Elongation	72	60 sec	
5. Amplification	72	7 min	1
6. Finish	4	∞	n/a

Primers:tTA PCR

Electrophoresis Protocol:

Name	Nucleotide Sequence (5' - 3')	Argarose: 2% in v:110		
1. tTA For	GCTGCTTAATGAGGTCGG	Estimated Running:Time: 60 min.		
2. tTA Rev	CTCTGCACCTTGGTGATC	Primer Combination	Band	Genotype
3.		tTA For + Rev	514 bp	Tet O (tTA)
4.			bp	
5.			bp	

Please size gel images to fit in this space

GENOTYPING BY PCR PROTOCOL
MUTANT MOUSE REGIONAL RESOURCE CENTER

sacoord@mmrc.org

800-910-2291 North America, +1-530-757-5710 International

Protocol / Gel Comments:

We digest mouse tails with Proteinase K overnight at 55 centigrade in a volume of 40uL, inactivate the Prot. K at 95 centigrade for 10 min. and then add 10mM TRIS, pH 8.0 to get a final volume of 400uL. We use 5uL of this crude genomic DNA for all our genotyping.

We run a separate 25uL genomic PCR reaction for huPMCA4b transgenic locus and another separate 25uL genomic PCR reaction for the Tet O (Tetracycline sensitive TransActivator transgenic locus. We use positive control tail DNA from previously genotyped mouse tail DNA as well as a 100bp-1,000 bp DNA ladder where the 500bp band is more intense than the rest of the ladder bands to assign genotypes.

Gel pictures: These are the results of our latest genotyping (2018-05-25).

