

**GENOTYPING BY PCR PROTOCOL
MUTANT MOUSE REGIONAL RESOURCE CENTER**

sacoord@mmrrc.org

800-910-2291 North America, +1-530-757-5710 International

Please provide the following information required for genetic analysis of your mutant mice.

Note to MAC users: to ensure your graphic can be viewed on a PC please follow the steps below when inserting the graphic into this document. DO NOT drag and drop or copy/paste the graphic into this document.

- Open the original graphic in the program that created it
- Choose File, Save As
- Select No Compression in the save options.
- Save as JPG or PNG or similar format that's compatible with both PC and Mac Word versions.
- Switch to Word, choose Insert, Picture, From File and choose the newly saved picture.

These instructions are very generic. The menu options for your graphics program may be different.

Donating Investigator/PI		
Yusuf Hannun		
Email		
Yusuf.Hannun@stonybrookmedicine.edu		
Institution		
Stony Brook University		
Address		
100 Nicolls Road		
City	State	Zip
Stony Brook	NY	11794-8155
Lab Contact		
Ashley Snider		
Email		
Ashley.Snider@stonybrookmedicine.edu		
Telephone	FAX	
631-444-7488	631-444-6313	
Strain Name		MMRRC Stock Number
P38 delta null		43950

GENOTYPING BY PCR PROTOCOL
MUTANT MOUSE REGIONAL RESOURCE CENTER

sacoord@mmrrc.org

800-910-2291 North America, +1-530-757-5710 International

NAME OF PCR: P38 Delta Null **MMRRC:** 0-CTR

Protocol: *(PCR protocol provided by Donating Investigator)*

Reagent/Constituent	Volume (μL)
Water	7.5
10x Buffer WE USE PROMEGA MASTER MIX (see comments)	12.5
MgCl ₂ (stock concentration is mM)	0
Betaine (stock concentration is 5M) <i>Optional</i>	0
dNTPs (stock concentration is 10mM)	0
DMSO <i>Optional</i>	0
Primer 1. (stock concentration is 10μM) P38_F Forward	1.0
Primer 2. (stock concentration is 10μM) P38_R Reverse	2.0
Primer 3. (stock concentration is 10μM) P38_NEO	1.0
Taq Polymerase 5Units/μL	0
DNA (50-200ng/ μL) extracted w/ "Qiagen DNeasy columns or other similar silica based kits"	1.0
<i>The total volume is auto-calculated based on volumes entered, right click the total and update field to show/recalculate the total volume.</i>	TOTAL VOLUME OF REACTION:
	25.000 μL

Comments on protocol:

•

Strategy:

Steps	Temp (°C)	Time (m:ss)	# of Cycles
1. Initiation/Melting HOT START? ☐	95	2.:00	1
2. Denaturation	95	:30	
3. Annealing steps 2-3-4 cycle in sequence	60	:30	40x
4. Elongation	72	1:00	
5. Amplification	72	10:00	1
6. Finish	4	∞	n/a

Primers:

Electrophoresis Protocol:

Name	Nucleotide Sequence (5' - 3')	Argarose: 1% V:
1. P38_F Forward	CCCTTGAGCCATAGATCCTGG	Estimated Running:Time: 30 min.
2. P38_R Reverse	CATGAGCTTGAGATGCTCTCT	Primer Combination
3. P38_NEO	GGCGATGCCTGCTTGCCGAAT	Band
4.		Genotype
5.		

Please size gel images to fit in this space

**GENOTYPING BY PCR PROTOCOL
MUTANT MOUSE REGIONAL RESOURCE CENTER**

sacoord@mmrrc.org


800-910-2291 North America, +1-530-757-5710 International

Protocol / Gel Comments:

For the PCR reaction we use GoTaq® Colorless Master Mix from Promega (catalog # M7133).

Gel pictures:

W H K W K K W _____ W H



Final lane is heterozygous positive control

W = wildtype

H = heterozygote

K = p38 delta null