

GENOTYPING PROTOCOL

MUTANT MOUSE RESOURCE & RESEARCH CENTER: UC DAVIS

mmrrc@ucdavis.edu

530-754-MMRRC

Protocol Name: CR1744 Trappc2l EXDEL

Protocol: GoTaq® G2 Colorless Master Mix(Promega)

| Reagent/Constituent | Volume (µL) |
|---|-----------------|
| Water | 4.5 |
| GoTaq® G2 Colorless Master Mix,2X | 7.5 |
| Primer 1. (stock concentration is 20µM) comF | 0.5 |
| Primer 2. (stock concentration is 20µM) wtR | 0.5 |
| Primer 3. (stock concentration is 20µM) mutR | 0.5 |
| DNA (example) extracted w/ "Qiagen DNeasy columns or other similar silica based kits" | 1.5 |
| TOTAL VOLUME OF REACTION: | 15.00 µL |

Comments on protocol:

- Protocol may work with other DNA extraction methods.

Strategy:

| Steps | Temp (°C) | Time (m:ss) | # of Cycles |
|--|-----------------|---------------------|-------------|
| 1. Initiation/Melting HOT START? <input type="checkbox"/> | 94 | 2:00 | 1x |
| 2. Denaturation | 94 | 0:10 | |
| 3. Annealing steps 2-3-4 cycle in sequence | 65 (↓1°C/cycle) | 0:30 | 10x |
| 4. Elongation | 68 | 2:00 | |
| 5. Denaturation | 94 | 0:15 | |
| 6. Annealing steps 5-6-7 cycle in sequence | 55 | 0:30 | 25x |
| 7. Elongation | 68 | 2:00 (↑20sec/cycle) | |
| 8. Finish | 4 | ∞ | n/a |

Primers:

Electrophoresis Protocol:

| Name | Nucleotide Sequence (5' - 3') | Agarose: 1.5% | V: 90 |
|---------------------|-------------------------------|---------------------------------|------------------|
| 1. CR-Trappc2l-comF | GGCTTAATGCTGCGACTCAAGC | Estimated Running Time: 90 min. | |
| 2. CR-Trappc2l-wtR* | GCTCACTCTCTGTCTGGAGTGCTTC | Primer Combination | Band (bp) |
| 3. CR-Trappc2l-mutR | GGAGTCAGAGGAAGTGCCAGGCT | 1 & 2, 1 & 3 | 351, 1037 |
| | | 1 & 3 | 355 |
| | | | wildtype |
| | | | mutant |

Allele Description: Exon 2 [ENSMUSE00000295747](#) and flanking splicing regions were constitutively deleted from the Trappc2l gene [ENSMUSG00000015013](#) using CRISPR Cas9 gene editing technology in mouse zygotes. Subsequent founders were backcrossed to C57BL6/N to produce sequence confirmed heterozygous animals.

*Run in Simplex

*wtR primer untested (ePCR verified)

